



- Class B
- LCD panel
- Pre & post vacuum
- Prion program
- Leakage, Helix, B.D. test program
- Customization program
- Data recorder (SD card) and printer
- Tank can be filled manually or piped-in
- Overheat protection
- Process evaluation system
- Pressure auto door lock
- Emergency stop button
- Pressure overload protection
- Chamber size 260X450
- Chamber volume : 24 L
- Wattage - 2770W

Modern Surgical Company

101/C, C. R. Avenue
Kolkata - 700 073
E-mail : modernsurgical84@gmail.com
Mob. : +91 98303 02017 / 93082 41770

www.modernsurgical.com