



Product Specification :

MODEL NO.	SA-232X
Chamber Size	230mmø x 410mm(D) Stainless Steel #304
Chamber Capacity	16 Liter
Overall	501(W) x 406(H) x 537(D)mm
Power consumption	220V - 240V / 1400W
Using Pressure and Temperature	UNWRAPPED / WRAPPED / LIQUIDS 118°C ~ 134°C(0.9~2.1 KG/CM ²)
Design Temperature	140°C
Safety Approval	CE-EMC, LVD Approved CE-MDD Approved Chamber design by ASME/U.S.A Standard
Standard Accessories	Stainless Steel #304 Sterilization Box./1pc.